

2010		Thermo-ARL Maxi Torch ICP's ( 3560, 3580, etc )
EPOND N°	OEM P/No.	Description
2050-0800N	915375	K-Type Concentric <b>Glass</b> Nebulizer, with <i>NebLink</i> sample tube
2050-M030-C1	TR30-C1	' <b>Meinhard</b> ' Concentric glass nebulizer, Type C
2051-0500	173259	Ryton Babbington Nebulizer
2050-0700NQ		Standard Concentric <b>Quartz</b> Nebulizer, with ' <i>NebLink</i> '
2050-0710N		Standard Concentric <b>Glass</b> Nebulizer, with ' <i>NebLink</i> '
2050-0710NT		High-Performance ' <b>Tyfoon</b> ' Concentric Glass Nebulizer, with ' <i>NebLink</i> '
2050-4254		<b>Boreal</b> 'NB-Class' ( <b>No-Break</b> ) Concentric Polymer Nebulizer, medium uptake
2550-2028		<b>PolyPro-ST</b> HF-res concentric nebulizer, with detachable sample uptake tube
2550-2024		<b>PFA-ST</b> HF-res concentric nebulizer, with detachable sample uptake tube
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		<b>Teflon</b> <i>MiraMist</i> HF-resistant Parallel-Path Nebulizer
3150-0720		<b>PEEK</b> <i>MiraMist</i> HF-resistant Parallel-Path Nebulizer
2010-0120	173142-0000	Maxi-Torch Conical Glass Spray chamber, for 6-mm nebs
2010-0140	173142-0001	Maxi-Torch MDSN Spray Chamber
2010-0160	139008-0000	Maxi Torch Conical Spray Chamber W/Impact Bead
2015-0175		' <b>Althea</b> ' Cyclonic Glass Spray Chamber, 12/5 socket, 4mm pumped drain
2015-0176		' <b>Libra</b> ' Cylonic Baffled Spray Chamber 12/5 Socket 4mm Drain
2045-0136		' <b>Electron</b> ' Mini-cyclonic spray chamber, volume 20 ml., RS13/5 outlet, 4mm drain
2045-0137		' <b>Positron</b> ' Baffled Mini-cyclonic spray chamber, volume 20 ml., RS13/5 outlet, 4mm drain
2011-0550		Cyclonic Glass Spray Chamber, 10mm Drain, 12/5 Outlet
2011-2929		' <b>Erebus</b> ' HF-resistant High-stability Cyclonic Spray Chamber
3180-3001		<b>PFA HF-resistant</b> Cyclonic Spray Chamber, RS12/5 Outlet
2010-0050	139009-0003	ARL Quartz 1-piece Maxi Torch 1.8mm injector
2010-050W		Wear Metals Maxi Torch
2010-0060		Maxi Torch with 20mm Shield
2010-0080	137000-0000	ICPQ Maxi Torch
2010-0090	173127-0000	Maxi Torch, Nitrogen, 3"Port 28/18 Ball
2010-0210		Maxi Torch Extension Tube
2010-0220		Maxi Torch Extension Tube, with Optical Port
2010-0250		Maxi Torch Bonnet
2010-0270	173251-0000	Maxi Torch Bonnet, Induction Coil M/S
2012-0400	173131	HF/CR Torch Body, for HF-resistant version of Maxi Torch
2012-0420		Alumina Injector Tip, for 2012-0400 Maxi Torch
2012-0440		Teflon Insert for 2012-0400 HF-res Maxi Torch
2012-0450		Adjustable Teflon Insert, for Maxi Torch
2012-0510		Teflon Insert for Maxi Torch MDSN Spray Chamber
EPOND No.		Accessory Description
QuDIN-1011		Mermet-Todoli QDIN adaptor, with Peri-Pump tubes and connectors
APEX-E-ARL		<b>APEX</b> Desolvator, with PolyPro-ST concentric nebulizer and a ' <i>Pluto</i> ' flow meter
2019-4260		' <b>XENA</b> ' micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
8100-2051		' <b>PEARL</b> ' Cleaning Device for nebulizers - <b>suits all types concentric glass nebs</b>
6808-2280		' <b>ISOBAR</b> ' Internal Standard Kit, HF-resistant, with tubing and connectors
2045-0700		' <b>SOLVIT</b> ' Argon Humidifier, with tubing and connectors

EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		<i>Cone-terminated</i> Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		<b>EPOND</b> 'RubeTubes', sample tube connectors, packet of 40
		Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a>
		<a href="http://www.e-pond.biz">www.e-pond.biz</a>