

2012		Thermo – ARL – Accuris, 3410, VGE Horizon and Mini-torch ICP's
EPOND N°	OEM P/No.	EPOND S.A. Parts List ... Description
2050-0820N	915375	K-Type ' <b>Caliber</b> ' 2ml/min. Concentric Nebulizer, <i>with NebLink</i> sample tube
2050-0770N	917319	A-Type 3ml/min. Concentric Nebulizer, <i>with NebLink</i> sample tube
2051-0100	915817	V-Groove MDSN Glass Nebulizer
2051-0500	173259	Ryton Babbington Nebulizer
2050-M030-C1	TR-30-C1	<b>'Meinhard'</b> Concentric glass nebulizer
2050-4250		<b>Boreal 'No-Break' HF-resistant</b> Concentric Polymer Nebulizer
2550-2028		<b>PolyPro-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		<b>PFA-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Burgener Teflon <b>MiraMist</b> Parallel-Path Nebulizer
3150-0720		Burgener PEEK <b>MiraMist</b> Parallel-Path Nebulizer
3150-0730		<b>Teflon T2100</b> 'Work-horse' Parallel-Path Nebulizer
2011-0330	173142-0002	ARL Mini Conical Glass Spray Chamber, 6-mm Nebs
2011-0330Q		Mini Torch Quartz Conical Spray Chamber, for 6mm nebs
2011-0340	173142-0003	Mini-Torch MDSN Conical Glass Spray Chamber
2011-0360	173142-0004	Mini-Torch Jacketed Conical Glass Spray Chamber
2011-2929		HF-resistant ' <b>Erebus</b> ' Cyclonic Spray Chamber, low RSD's
2051-026T		PTFE HF-resistant Cyclonic Spray Chamber
2011-0550		<b>'Althea'</b> Cyclonic Spray Chamber, 10mm Drain, 12/5 Outlet
2011-055J		Jacketed <b>'Althea'</b> Cyclonic Spray Chamber, 10mm Drain, 12/5 Outlet
2011-0175Q		Jacketed Quartz <b>'Althea'</b> Cyclonic Spray Chamber 10mm Drain 12/5 Socket
2015-0175		<b>'Althea'</b> Cyclonic Glass Spray Chamber, 12/5 socket, 4mm pumped drain
2015-0176		<b>'Libra'</b> Cylonic Baffled Spray Chamber 12/5 Socket 4mm Drain
2045-0136		<b>'Electron'</b> Mini-cyclonic spray chamber, volume 20 ml., RS13/5 outlet, 4mm drain
2045-0137		<b>'Positron'</b> Baffled Mini-cyclonic spray chamber, vol. 20 ml., RS13/5 outlet, 4mm drain
2011-0300	173463-0000	Standard ARL Quartz Mini Torch, Inj. ID=1.5mm, 4mm arms
2011-030U		ARL Quartz Mini Torch, Inj. ID=1.0mm
2011-6831		Semi-Demountable HF-resistant Mini-Torch, with alumina injector
2011-0310	173473-0000	ARL Mini Torch HF/CR Torch Body, L = 90mm., fits SDT No. 2011-6831
2011-0320	170692-0000	ARL Mini Torch, for Nitrogen apps
2011-0350	173252-0000	Quartz Bonnet, for ARL Mini Torch
2011-0420	144860	Mini Torch - Accuris, Torch Mount
2809-0050		HF-res. PP Connector, Mini-Torch to Spray Chamber RS13/5
2011-0480	144857	HF-resistant Mini Torch Alumina Injector, fits SDT No. 2011-6831
2011-0500		Mini Torch, High Solids version
2011-0540	170693	Mini Torch Argon Snout N2
<b>EPOND No.</b>		<b>Accessory Description</b>
2019-4260		<b>XENA'</b> micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
8100-2051		<b>'PEARL'</b> Cleaning Device for glass nebulizers - <b>suits all types concentric glass nebs</b>
6808-2280		<b>'ISOBAR'</b> Internal Standard Kit, HF-resistant, with tubing and connectors
2045-0700		<b>'SOLVIT'</b> Argon Humidifier, with tubing and connectors
<b>EPOND N°</b>	<b>OEM P/No.</b>	<b>Description, tubes</b>
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
ERT-040		<b>EPOND 'RubeTubes'</b> , sample tube connectors, packet of 40
DLK-075		<i>DrainLink</i> nebulizer adaptor and tube with 0,75mm ID, pkt 5
NLK-030		<i>Neb-Link</i> sample tubes, packet of 10, ID=0,30 mm, L=70cm.
ICP-standards		<i>Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a></i>
		<i>... custom standards supplied - according to your specifications ...</i>
		 
		<a href="http://www.e-pond.biz">www.e-pond.biz</a>