

2010		Components for LECO ICP-MS
EPOND No.	OEM P/No.	Description
2050-0710N		Standard Concentric Nebulizer, with Neblink sample tube
2050-0770Q		A-Type 3ml/m Quartz Concentric Nebulizer
2050-0710NT		High-Performance 'Tyfoon' Concentric Glass Nebulizer
2050-M030-C1	TR30-C1	'Meinhard' Concentric glass nebulizer, Type C
2050-4151		'Zefyr' HF-resistant Concentric Polymer Nebulizer, low-uptake version
2050-4251		'Boreal' HF-resistant Concentric Polymer Nebulizer, medium-uptake version
2550-2000		HF-resistant MicroFlow PFA Concentric Nebulizer, 50 ul/min uptake version
2550-2002		HF-resistant MicroFlow PFA Concentric Nebulizer, 100 ul/min uptake version
2550-2003		HF-resistant MicroFlow PFA Concentric Nebulizer, 200 ul/min uptake version
2550-2005		HF-resistant MicroFlow PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		<b>PolyPro-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		<b>PFA-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2045		Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min.
2550-2043		Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min.
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
		<b>Capillaries with integrated probes - available for all common autosamplers</b>
3150-0710		<b>Burgener</b> PTFE <b>MiraMist</b> Parallel-Path Nebulizer
3150-0720		<b>Burgener</b> PEEK <b>MiraMist</b> Parallel-Path Nebulizer
2015-7175		Standard ' <b>Althea</b> ' 50ml. Cyclonic glass spray chamber
2015-7176		Standard ' <b>Libra</b> ' 50ml. Baffled Cyclonic glass spray chamber
2045-7136		Standard ' <b>Electron</b> ' 22ml. Mini-cyclonic glass spray chamber
2045-7137		Standard ' <b>Positron</b> ' baffled Mini-cyclonic glass spray chamber
2045-7929		' <b>Erebus</b> ' HF-resistant cyclonic spray chamber, with RS13/5 outlet
EP-4157-1000	ES-4157-1000	ESI ' <b>PC3</b> ' Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply
EP-9100-5400	ES-9100-5400	' <b>PC3-LT</b> ' Peltier-cooled Quartz Spray Chamber, low - temp. version ( -20 deg. C. )
2090-0380-Q		Quartz Extender Arm RS13/5
2090-0050	701-030	1-piece Quartz Torch with 2.5mm Injector ID
2090-0070	700-670	1-piece Quartz Torch with 1.5mm Injector ID
2090-0090		Semi Demountable Quartz Torch, with alumina injector for samples with HF
2090-0092		Semi Demountable Quartz Torch Body, with 6mm side arms
2066-0513		Semi Demountable Quartz Torch Body, side arms 6mm and 4mm
2066-0777		Semi Demountable Torch, 6mm and 4mm side arms, with Quartz Injector
<b>EPOND N°</b>		<b>Other Accessories</b>
APEX-Q-LECO		<b>APEX-Q</b> Desolvator, with 2 PFA concentric nebulizers + 'Sedna' Flow Meter
2019-4260		' <b>XENA</b> ' micro-flow meter / aspiration monitor, range : 5 - 2000 ul/min., mains powered
ACM-E2		ESI Nafion <b>Active Membrane Device</b> , with DC Power Supply
EP-4599-3002	ES-4599-300A	<b>Spiro TMD</b> Membrane Desolvation Module ( latest model of ACM ), with Power Supply
6808-2280		' <b>ISOBAR</b> ' Internal Standard Kit, HF-resistant, with tubing and connectors
8100-2051		' <b>PEARL</b> ' Cleaning Device for concentric glass nebulizers - <b>suits all types</b>
<b>EPOND N°</b>	<b>Leco No.</b>	<b>Description of Cones</b>
LE1001-Ni	703-633	Nickel Sampler Cone

LE1002-Ni	703-634	Nickel Skimmer Cone
LE1001-Ni/Cu	703-633	Nickel Sampler with copper base
LE1002-SS	710-804	Ni Skimmer Cone, electro-polished
LE1006-Pt	703-270	Platinum Sampler Cone
LE1008-Pt		Platinum Skimmer Cone
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe
		<i>Note : <b>Detector EPM-14833-H</b> available, price CHF 2795,- / Euros 1870,-</i>
<b>EPOND N°</b>	<b>OEM P/No.</b>	<b>Description, tubes</b>
EPP-040		IDEX High-Precision 4-Channel Peristaltic Pump
<i>See web site</i>		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
<i>See web site</i>		<i>Cone-terminated</i> Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		<b>EPOND</b> 'RubeTubes', sample tube connectors, packet of 40
		<i>Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a></i>
		<b>reference solutions made to measure - on request</b>
		<b><a href="http://www.e-pond.biz">www. e-pond. biz</a></b>