

2012		EPOND Parts list
EPOND N°	OEM P/No.	Description, accessories for Leeman ICP instruments
2050-0710N		C-Type standard Concentric Glass Nebulizer, 2 ml/min., with <i>NebLink</i> tube
2050-0710NT		' Tyfoon ' Hi-performance Concentric Nebulizer, with <i>NebLink</i> sample tube
2050-M030-C1	TR30-C1	' Meinhard ' Concentric glass nebulizer, Type C
2050-4253		'NB-Class' (No-Break) Boreal Concentric Polymer Nebulizer, medium uptake version
2050-4353		'NB-Class' (No-Break) Capital Concentric Polymer Nebulizer, high uptake version
2550-2005		HF-resistant MicroFlow PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Teflon Burgener MiraMist Parallel-Path Nebulizer
3150-0720		PEEK Burgener MiraMist Parallel-Path Nebulizer
2055-0550		Scott-type Coaxial standard glass spray chamber
2052-0300		Adaptor, PTFE 35/6-mm for Scott chamber – glass nebulizer
2055-0570		'Althea' Cyclonic glass spray chamber, for standard nebulizers with 6-mm diameter
2055-0580		Cyclonic spray chamber, for Hildebrand Grid nebulizer
2055-2929		HF-resistant 'Erebus' Cyclonic spray chamber
2055-0640	318-00064	Aqueous Axial DRE quartz torch
2055-0650	318-00065	Aqueous Dual-View standard DRE 1-piece quartz torch
2055-0670	318-00067	Purge Tube, Dual-View DRE
2055-0680	318-00068	Radial DRE Purge Tube
2055-0690	318-00069	Axial-View Demountable Torch. Aqueous
2055-0760	318-00076	Axial-View Torch, Aqueous with injector ID 1,5mm
EPOND N°		Flow Meters, Desolvators and Other Accessories
APEX-Q-LE		APEX Desolvator, with PolyPro concentric nebulizer and XENA uptake meter
2019-4260		' XENA ' micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply
8100-2051		' PEARL ' Cleaning Device for concentric glass nebulizers – <i>suits all types !</i>
2045-0700		' SOLVIT ' Argon Humidifier, with tubing and connectors - <i>prevents salting up of nebulizers</i>
6808-2280		' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
ERT-040		EPOND ' <i>RubeTubes</i> ', sample tube connectors, packet of 40
DLK-075		<i>DrainLink</i> nebulizer adaptor and tube with 0,75mm ID, pkt 5
NLK-030		<i>Neb-Link</i> sample tubes, packet of 10, ID=0,30 mm, L=60cm.
		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards
		... custom standards supplied - according to your specifications ...
		www.e-pond.biz