

2010		EPOND S.A. Price List, PE Optima Radial View ICP
EPOND N°	OEM P/No.	Description
2050-0700N	N810-2011	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-0700NQ		Standard Concentric Quartz Nebulizer, with ' <i>NebLink</i> ' sample tube
2050-0710N	0047-2022	Hi-flow Concentric Glass Nebulizer, with <i>NebLink</i>
2050-0710N	N068-1574	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-M030-C1	TR30-C1	' Meinhard ' Concentric glass nebulizer, Type C
2050-0710NT		High-Performance ' Tyfoon ' Concentric Glass Nebulizer, with <i>NebLink</i>
2050-4250		Boreal 'No-Break' HF-resistant Concentric Polymer Nebulizer
2550-2005		HF-resistant MicroFlow PFA Concentric Nebulizer, 360 ul/min uptake version
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Burgener Teflon MiraMist Parallel-Path Nebulizer
3150-0720		Burgener PEEK MiraMist Parallel-Path Nebulizer
3150-0730		Teflon T2100 'Work-horse' Parallel-Path Nebulizer
2052-0300		PTFE Endcap for CG Nebulizer, 35/6 mm adaptor
2030-0150		Optima 3000 Radial Spray Chamber, Top Port Angled
2030-0160		Scott (Coaxial) Spray Chamber, Optima 3000 Radial
2030-0190		Cyclonic Baffled Spray Chamber for Optima 3000 radial
2030-0200		MLN Cyclonic Spray Chamber, Optima 3000 Radial
2030-0194		' Libra ' baffled cyclonic glass spray chamber, 50 ml. With 11-mm ID output tube
6609-0222		' Libra ' baffled cyclonic glass spray chamber, 50 ml. With RS13/5 output
2031-0350	6808-0204	Spray-Chamber Torch adaptor for cyclonic spray chambers
EP-4157-1000		PC3 Peltier-cooled Cyclonic Spray Chamber, with 12-Volt Power Supply
EP-9100-5700		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2031-0040	N069-0568	Optima 3000 Radial Quartz Torch
2031-004W		Optima 3000 Radial Torch, for Wear Metals application
2031-0050	N069-5295	Torch Bonnet for Optima 3000 radial, before sept. 1994
2031-0130	N069-5456	Optima 3000 Radial Qtz Torch Bonnet - post Sept. 1994
2032-0330	N068-1631	Qtz Sample Injector Optima 3000 Radial, ID=1.2mm
2032-0430	N069-5448	Qtz Sample Injector Optima 3000 Radial, ID=1.6mm
2032-0400	N069-5441	Qtz Sample Injector Optima 3000 Radial, ID=2.0mm
2032-0510	N058-2186	Alumina Injector, ID=0.85mm
2032-0520	N058-2185	Alumina Injector, ID=1.5mm
2032-0530	N058-2184	Alumina Injector, ID=2.0mm
EPOND P/N°		Other Sample Introduction Accessories
APEX-E-PE		APEX-E Desolvator, with PolyPro-ST nebulizer and PLUTO uptake meter
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module with Power Supply
2019-4260		' XENA ' micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
8100-2051		' PEARL ' Cleaning Device for concentric glass nebulizers
6808-2280		' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors
2045-0700		' SOLVIT ' Argon Humidifier, with tubing and connectors
EPOND N°	OEM P/No.	Description, tubes
EPP-040		IDEX High-Precision 4-Channel Peristaltic Pump
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40
EST-030		PTFE Sample Tube, OD=1,3mm ID=0,30mm L=10m.
EST-050		PTFE Sample Tube, OD=1,3mm ID=0,50mm L=10m.
DLK-075		<i>DrainLink</i> nebulizer adaptor and tube with 0,75mm ID, pkt 5
ICP-standards		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards

		<i>... custom standards supplied - according to your specifications ...</i>
		<i>www. e-pond. biz</i>