

2012		EPOND Parts List, PE Plasma 400 /1000/2000
EPOND N°	OEM P/No.	Descriptions
2050-0700N	N810-2011	Standard ' Caliber ' Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-M030-C1	TR30-C1	'Meinhard' Concentric glass nebulizer, Type C
2050-0700NQ		Standard Concentric Quartz Nebulizer, with ' <i>NebLink</i> ' sample tube
2050-0720N	0047-2022	Hi-flow Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-0710N	N068-1574	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-0710NT		High-Performance ' Tyfoon ' Concentric Glass Nebulizer, with <i>NebLink</i>
2050-4250		Boreal 'No-Break' HF-resistant Concentric Polymer Nebulizer
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Burgener Teflon <i>MiraMist</i> Parallel-Path Nebulizer
3150-0720		Burgener PEEK <i>MiraMist</i> Parallel-Path Nebulizer
2050-0920		HF-resistant Fixed Crossflow Nebulizer, for Scott type chamber

3304-5791P		Jacketed ' Althea ' Cyclonic Glass Spray Chamber for volatile samples
2030-0150P		Glass Scott Spray Chamber, Top Port Angled
2030-0160P		Glass Scott (Coaxial) Spray Chamber
2030-0190P		Libra' Cyclonic Baffled Spray Chamber
ES-4157-1000		PC3 Peltier-cooled Quartz Spray Chamber, with Power Supply

2031-0030	N058-0528	Plasma 400 / 1000 / 2000 Quartz Torch-Type II, 3 slots
2031-003W	N068-1690	Plasma 400 / 1000 / 2000 Torch Type II, for Wear Metals
2031-0090	N058-0376	PE Plasma 40 Demountable Torch Tube, 2-slots
2031-0240	N058-2235	PE Plasma 400/1000/2000 Torch Quartz Bonnet
2031-0350	N058-2248	Plasma 400 / 1000 / 2000 Injector Support Adaptor with O-Rings
2032-0300	N068-1631	PE 400 Quartz Injector Tube ID=1.2mm
2032-0310	0047-3292	PE 400 Quartz Injector Tube ID=1.6mm
2032-0250		Quartz Injector 2.0 x 109
2032-0300		PE 6500 Quartz Sample Injector Tube ID=1.2mm, Type I
2032-0310	N058-0449	PE 6500 Quartz Sample Injector Tube ID=1.6mm, Type I
2032-0320		PE 6500 Qtz. Sample Injector Tube ID=2.0mm, Type I
2032-0600	N068-1234	PE Plasma 40/400 Purge Extension
2032-0660	N058-1497	Plasma 1000/2000 Purge Extension Window, UV

EPOND P/N°		Sample Introduction Accessory description
SPIRO-TMD	EP-4599-3002	'Spiro-TMD' autonomous Membrane Separator, for aerosol solvent removal
APEX-Q		'APEX-Q' desolvator , for ~ 10x sensitivity improvement

2019-4260		'XENA' micro-flow meter / aspiration monitor, range : 5 – 1999 ul/min.
QuDIN-1011PE		Mermet-Todoli QDIN adaptor, with Peri-Pump tubes and connectors
2045-0700		'SOLVIT' Argon Humidifier, with tubing and connectors
8100-2051		'PEARL' Cleaning Device for concentric glass nebulizers – <i>suits all types</i>
6808-2280		'ISOBAR' Internal Standard Kit, HF-resistant, with tubing and connectors
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40
DLK-075		DrainLink nebulizer adaptor and tube with 0,75mm ID, pkt 5
ICP-standards		<i>Certified ICP / AAS Standard solutions - www.e-pond.biz/standards</i>
		<i>www.e-pond.biz</i>