

2010		EPOND Parts List for the Spectro CIROS, GENESIS and ARCOS
EPOND P/N°	OEM P/No.	Component descriptions - radial and axial
EP-2050-0700N	76060510	C-Type Concentric Glass Nebulizer, with 'NebLink' sample tube
EP-2050-0700NQ		Standard Concentric <b>Quartz</b> Nebulizer, with 'NebLink' sample tube
EP-2050-0710NT		High-Performance ' <b>Tyfoon</b> ' Concentric Glass Nebulizer, with sample tube
2050-M030-C1	TR30-C1	<b>'Meinhard'</b> Concentric glass nebulizer, Type C
EP-2050-4251		<b>'Boreal'</b> Medium-Uptake ' <b>No-Break</b> ' concentric polymer nebulizer
EP-2050-4351		<b>'Capital'</b> High-Uptake ' <b>No-Break</b> ' concentric polymer nebulizer
EP-2550-2002		HF-resistant PFA MicroFlow Concentric Nebulizer, 100 ul/min uptake version
EP-2550-2005		HF-resistant PFA MicroFlow Concentric Nebulizer, 400 ul/min uptake version
EP-2550-2028		<b>PolyPro-ST</b> HF-res concentric nebulizer, with detachable sample uptake tube
EP-2550-2024		<b>PFA-ST</b> HF-res concentric nebulizer, with detachable sample uptake tube
EP-2030-2800	ES-2030-28	<b>PFA-HS</b> High-Solids HF-res concentric nebulizer, with detachable sample uptake tube
EP-2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
EP-2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
EP-5017-0750-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,8 mm
EP-5017-0500-100		CETAC ASX-500 Autosampler Probe with Capillary for PFA-ST nebulizer, ID = 0,5 mm
EP-3150-0710		Teflon <b>Burgener MiraMist</b> Parallel-Path Nebulizer
EP-3150-0720		PEEK <b>Burgener MiraMist</b> Parallel-Path Nebulizer
3150-0730		Teflon <b>Burgener T2100</b> 'Work-horse' Parallel-Path Nebulizer
EP-2050-0925	75060502	Fixed Cross-Flow for Spectro, 34mm Inner Diameter
EP-2051-0205	48205036	Spectro Modified Lichte Nebulizer
EP-2015-0175		<b>'Althea'</b> Cyclonic Glass Spray Chamber, 12/5 socket, 4mm pumped drain
EP-2027-5161		<b>'Electron'</b> Mini-Cyclonic Glass spray chamber, vol. 20 ml.
EP-2027-5166		<b>'Positron'</b> Baffled Glass Mini-Cyclonic spray chamber, vol. 20 ml.
EP-2075-0120	48105003	Spectro Small Conical Spray Chamber, with Impact Sphere
EP-2075-0130	48105051	Spectro Glass Spray Chamber, Double Pass 90° Exit
EP-2075-0130	4810501	Spectro Quartz Spray Chamber, Double Pass 90° Exit
EP-3167-1001	ES-3167-1001	Quartz cyclonic spray chamber, can be used with PC3 Peltier cooling unit
EP-2075-0136		<b>'Erebus'</b> HF-resistant Cyclonic spray chamber
EP-2075-0140	48105006	Spectro Large Double Pass Spray Chamber, RS12/5 outlet
EP-2100-2477		HF-resistant PFA spray chamber kit with PP adaptor for Spectro X-Flow nebulizer
EP-2256-0470		47-mm diameter PFA spray chamber for ICP-OES, with 'PureCap' end-cap
EP-2256-5470		47-mm PFA spray chamber for Spectro ICP-MS, for use with ESI torch
EP-2100-2478		PolyPropylene adaptor for use of Spectro X-Flow nebulizer with 47-mm sp. Chamber
EP-4157-1000		<b>PC3</b> Peltier-cooled Quartz Spray Chamber, with 12-Volt DC Power Supply
EP-9100-5700		<b>PC3-LT</b> (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
EP-2250-0470	ES-2250-0470	PFA Sapphire HF-resistant Kit for Spectro CIROS, with PFA-ST nebulizer
EP-2075-0170	48105062	Spectro Elbow Adapter, EOP torch to Chamber
EP-2075-0180	48105002	Spectro Elbow Adapter Slurry, EOP Torch to Chamber
EP-2075-0300		Spectro Drain/Pressure Lock
EP-2075-0350		Spectro U-Drain Trap, 13mm OD
EP-2075-1702	ES-1702-0000	Quartz Torch for demountable axial-view CIROS
EP-2075-0210	48105008	Spectro Quartz Injector/Aerosol Tube 1.5 mm
EP-2075-0220	48105010	Spectro Quartz Injector/Aerosol Tube 1.8 mm
EP-2075-0230	48105009	Spectro Quartz Injector/Aerosol Tube 2.5
EP-1980-0180	ES-1980-0180	Sapphire Injector in PFA adaptor, ID = 1,8 mm

EP-1980-0300	ES-1980-0300	Sapphire Injector in PFA adaptor, ID = 3 mm
EP-1981-0200	ES-1981-0200	<b>Platinum</b> Injector in PFA adaptor, ID = 2 mm
EP-2075-0400	48105052	Spectro EOP Flared End Quartz Torch 2.5mm injector, 4mm side arms
EP-2075-0400F	75160526	Spectro EOP Flared End Quartz Torch 2.5mm injector, with fittings
EP-2075-0400LL		Spectro EOP Flared End Long-Life Quartz Torch
EP-2075-0410	48105052A	Spectro EOP Fixed Flared End Quartz Torch, 1.8mm ID injector
EP-2075-0420	48206007	Spectro EOP Demountable Torch Tube Assembly
EP-2075-0430	75160538	PTFE Housing for Demountable Torch
EP-2075-0440		Spectro EOP Demountable Flared End Torch
EP-2075-0450	75160545	Complete Demountable Torch
EP-2075-0460	48105058	Quartz Injector Tube, 2.5mm ID
EP-2075-0480	75160537	Spectro EOP Fixed Sheath Gas System Torch
EP-2075-0550		Gas Connectors with Part Nos. 44302072 and 443302077
EP-2075-0930		Spectro Replacement O-Ring for Endcap
EP-1902-0210		1-piece Quartz Torch for radial-view CIROS, with 1,0mm injector
EP-2075-1903		1-piece Quartz Torch for radial-view CIROS, with 1,8mm injector
EP-1902-0225		1-piece Quartz Torch for radial-view CIROS, with 2,5mm injector
EP-2075-1992		Complete Demountable Torch with quartz injector, for <i>radial-view</i> Spectro CIROS
EP-2075-1902	ES-1902-0000	Quartz Torch for demountable <i>radial-view</i> CIROS torch P/No. EP-2075-1992
EP-1983-2102		Quartz Injector with 1,0mm ID, for use with EP-2075-1992 Quartz radial Torch
EP-1983-2152		Quartz Injector with 1,5mm ID, for use with EP-2075-1992 Quartz radial Torch
EP-1983-2202		Quartz Injector with 2,0mm ID, for use with EP-2075-1992 Quartz radial Torch
EP-1980-00SA		Sapphire Injector for use with EP-2075-1992 Quartz Torch, ID range 1,0 - 2,0 mm.
EP-1981-00PT		<b>Platinum</b> Injector for use with EP-2075-1992 Quartz Torch, ID range 1,0 - 2,0 mm.
<b>EPOND No.</b>		<b>Desolvators and other accessories</b>
APEX-E-SP		<b>APEX-E</b> Desolvator, with PolyPro-ST nebulizer and 'PLUTO' Flow Meter
EP-4599-3002	ES-4599-300A	<b>Spiro TMD</b> Membrane Desolvation Module, with Power Supply
EP-2019-4260		'XENA' micro-flow meter / aspiration monitor, range : 5-1999 ul/min., AC power
EP-8100-2051		<b>'PEARL'</b> Cleaning Device for concentric glass nebulizers - <b>for all types !</b>
EP-6808-2280		<b>'ISOBAR'</b> Internal Standard Kit, HF-resistant, with tubing and connectors
EP-2045-0700		<b>'SOLVIT'</b> Argon Humidifier, with tubing and connectors
EP-5016-3250		Spectro AS400 Autosampler Probe, 0.25-mm ID, ~ 100 ul/min.
<b>EPOND No.</b>	<b>OEM P/No.</b>	<b>Descriptions of tubing accessories</b>
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		<b>EPOND 'RubeTubes'</b> , Tygon-core sample tube connectors, packet of 40
PRT-040		<b>EPOND 'RubeTubes'</b> , PharMed sample tube connectors for organics, packet of 40
DLK-075		DrainLink spray chamber drain adaptor and tube with 0,75mm ID, pkt 5
		<b>Autosamplers and ICP / AAS Standard and Reference Solutions</b>
SC-8000-XXXX		Complete range of ESI SC-series autosamplers - please ask for details
ICP-standards		Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a>
		<b>... custom standards supplied - according to your specifications ...</b>
		<b><a href="http://www.e-pond.biz">www.e-pond.biz</a></b>