

EPOND Part.No	OEM Part-No.	Description
<b>Nebulizers/Spray chambers/accessories</b>		
<b>Quartz Nebulizers</b>		
EP 2050-SILQ-50-05		MEINHARD® SilQ+ Ultra high purity quartz nebulizer
EP 2050-TQP50-A05		MEINHARD® TQ+ Ultra high purity quartz nebulizer, 0.5 ml/min uptake, 50 PSI
EP 2050-MQ50-A05		MEINHARD® Type A concentric quartz nebulizer, 0.5 ml/min uptake, 50 PSI
<b>Glass Nebulizers</b>		
EP 2050-TRP050-05	1161870	MEINHARD® TR+ glass nebulizer, 0.5 ml/min uptake
EP 2050-M050-A05		MEINHARD® Type A glass nebulizer, 0.5 ml/min uptake, high sensitivity
<b>PFA Nebulizers</b>		
ES-2040-73	1317090	MicroFlow PFA-ST nebulizer
ES-2050-52	1162030	PFA-LC, MicroFlow PFA nebulizer
<b>PFA Nebulizers with Integrated Probes</b>		
EP 2020-3503-080		PFA-20 MicroFlow nebulizer with integrated autosampler probe, 20µl/min uptake
EP 2000-3503-080		PFA-50 MicroFlow nebulizer with integrated autosampler probe, 50µl/min uptake
EP 2002-3503-100		PFA-100 MicroFlow nebulizer with integrated autosampler probe, 100µl/min uptake
EP 2003-3503-100		PFA-200 MicroFlow nebulizer with integrated autosampler probe, 200µl/min uptake
EP 2005-3503-100		PFA-400 MicroFlow nebulizer with integrated autosampler probe, 400µl/min uptake
EP 2020-3505-080		PFA-20 MicroFlow nebulizer with integrated autosampler probe, 20µl/min uptake
EP 2000-3505-080		PFA-50 MicroFlow nebulizer with integrated autosampler probe, 50µl/min uptake
EP 2002-3505-100		PFA-100 MicroFlow nebulizer with integrated autosampler probe, 100µl/min uptake
EP 2003-3505-100		PFA-200 MicroFlow nebulizer with integrated autosampler probe, 200µl/min uptake
EP 2005-3505-100		PFA-400 MicroFlow nebulizer with integrated autosampler probe, 400µl/min uptake
<b>PFA Sample Capillaries for ST Nebulizers</b>		
EP 2550-2045		Self-aspirating sample capillary for ST Nebulizer, 0.15mm ID(red marker), 20µL/min
EP 2550-2043		Self-aspirating sample capillary for ST Nebulizer, 0.20mm ID(purple marker), 50µL/min
EP 2550-2042		Self-aspirating sample capillary for ST Nebulizers, 0.25mm ID(green marker), 100µL/min
EP 2550-2047		Self-aspirating sample capillary for ST Nebulizers, 0.30mm ID(yellow marker), 200µL/min
EP 2550-2041		Self-aspirating sample capillary for St Nebulizers, 0.50mm ID(orange marker), 400µL/min
EP 2550-2044		Self-aspirating sample capillary for ST Nebulizers, 0.80mm ID(blue marker), 700µL/min
EP 2550-2049		Self-aspirating sample capillary for ST Nebulizers, 1.00mm ID(gray marker), 1mL/min
<b>Cyclonic Spray Chambers</b>		
EP 2048-0300	1317080	Quartz cyclonic spray chamber
EP 2048-0320	1320260	PFA cyclonic spray chamber
<b>Mini Cyclonic Spray Chamber</b>		
EP 2048-0400		Quartz mini-cyclonic spray chamber(20ml) with baffle
<b>Torches/shields</b>		
<b>Quartz Torch</b>		
EP 2048-0100	1230790	Quartz torch
<b>Torch Holder</b>		
EP 2048-1000	1231000	Torch Holder

**Injectors and accessories****Quartz Injectors**

EP 2048-0150		1.5mm quartz injector
EP 2048-0160	1305640	2mm quartz injector
EP 2048-0170	1305600	2.5mm quartz injector

**Sapphire Injectors**

EP 2048-0510		1.5mm sapphire injector
EP 2048-0500	1323600	1.8mm sapphire injector
EP 2048-0520		2.5mm sapphire injector

**Platinum Injectors**

EP 2048-0190		1mm platinum injector
EP 2048-0200		1.5mm platinum injector
EP 2048-0180	1321200	2mm platinum injector

**Cones****Sampling Cones**

VG1021-Ni	3600812	Ni sampler cone, Nicone TM, base for iCAP and VG ExCell.
EP-01289	3601289	Pt sampler cone(15.2mm) for iCAP Q and VG, all X

**Skimmer Cones**

EP-11870	1311870	Ni skimmer cone for iCAP Q, insert ready
EP-24540	1324540	Pt skimmer cone for iCAP Q, insert ready

**Skimmer Cone Inserts**

EP-1880	1311880	Nickel Skimmer Insert-2.8-High Sensitivity
EP-18470	1318470	Nickel Skimmer Insert-3.0-Standard for A and C-High Matrix
EP-18480	1318480	Nickel Skimmer Insert-3.5-High Matrix Tolerance

**Sampling Cone Gasket**

EP 2048-9000	1310900	Graphite Gasket for Sampler Cone with interface holes, 3 pack
EP 2048-4382	3004382	Graphite Gasket for Sampler Cone, 3 pack

**Cone Cleaning Kit&Detergent**

EPMNTKIT1		Cone Cleaning Kit
EPCitranox		Citranox Detergent

**Accessories****Gas Line**

EP 2502-1000		Additional gas line
--------------	--	---------------------

**Autosampler Probes**

ES-5037-3255-150		0.25mm ID SC-DX autosampler probe, carbon fiber support, 150cm
ES-5037-3505-150		0.5mm ID SC-DX autosampler probe, carbon fiber support, 150cm

**PFA Elbow**

ES 3150-0073-L	1320270	PFA elbow with add gas port
----------------	---------	-----------------------------

**Copper Load Coils**

EP 2048-5260		Copper load coil for Thermo(VG) iCAP Q
EP 2048-5260-AG		Silver plated copper load coil for Thermo(VG) iCAP Q