

2012		EPOD Parts List : Thermo Elemental X-series
EPOD N°	OEM P/No.	Descriptions
2050-0710N	1040894	Standard ' <b>Caliber</b> ' Concentric Glass Nebulizer, with ' <b>NebLink</b> ' sample tube
2050-0700NQ		Standard Concentric <b>Quartz</b> Nebulizer, with <b>PFA argon and sample connectors</b>
2050-070QDAC-FI	4600294-03	Concentric Glass Nebulizer, type C, 35psi, 1 ml/min.
2050-0710NT		High-Performance ' <b>Tyfoon</b> ' Concentric Glass Nebulizer, nominal uptake 1 ml/min.
2050-0704NT		High-Performance ' <b>Tyfoon</b> ' Concentric Glass Nebulizer, nominal uptake 0,4 ml/min.
2050-M030-C1	TR30-C1	' <b>Meinhard</b> ' Concentric glass nebulizer, Type C, with <b>PFA connectors</b>
2050-M030-A3		' <b>Meinhard</b> ' Concentric glass nebulizer, Type A3, with <b>PFA argon and sample line connectors</b>
2050-M030-K		' <b>Meinhard</b> ' Concentric glass nebulizer, Type K, with <b>PFA argon and sample line connectors</b>
2050-4251		'NB-Class' ( No-Break ) <b>Boreal</b> Concentric HF-res Nebulizer, medium uptake version
2050-4151		'NB-Class' ( No-Break ) <b>Zefyr</b> Concentric HF-res Nebulizer, low uptake version
2550-2000		HF-resistant <b>MicroFlow</b> PFA Concentric Nebulizer, 50 ul/min uptake version
2550-2002		HF-resistant <b>MicroFlow</b> PFA Concentric Nebulizer, 100 ul/min uptake version
2550-2003		HF-resistant <b>MicroFlow</b> PFA Concentric Nebulizer, 200 ul/min uptake version
2550-2005		HF-resistant <b>MicroFlow</b> PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		<b>PolyPro-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		<b>PFA-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2045		Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min.
2550-2043		Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min.
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
2046-0001		Obstruction Removal Kit for ST-type nebulizers, with fused silica capillaries
		<b>Capillaries with integrated probes - available for all autosamplers</b>
3150-0710		<b>Burgener</b> PTFE <b>MiraMist</b> Parallel-Path Nebulizer
3150-0720		<b>Burgener</b> PEEK <b>MiraMist</b> Parallel-Path Nebulizer
2045-0220	3204313	Jacketed Glass Spray Chamber with Impact Sphere
2045-0130		Conical Spray Chamber, with impact bead - older type
2045-0150	4600508	Conical Spray Chamber, with impact bead X series 2
2045-0810	3600170	Conical Spray Chamber, for Peltier system - standard for X-series
2045-0132		PTFE Drain Adaptor, for Conical spray chamber P/No. 2045-0130 ( old type )
2045-0134		PTFE Drain Adaptor, for Conical spray chamber P/No. 2045-0810
2045-0140		VG Horizon Conical Spray Chamber for 6-mm OD Nebulizers
ES-4190-60-29		<b>PC3</b> Sample Introduction Kit for XSeries2 with Quartz Spray Chamber and PFA-ST nebulizer
ES-9190-60-29		<b>PC3-LT</b> Sample Introduction Kit for XSeries2 with Quartz Spray Chamber and PFA-ST nebulizer
2045-1261		' <b>Althea</b> ' Cyclonic Spray Chamber, RS13/5 output
2045-1263		' <b>Libra</b> ' Baffled Cyclonic Spray Chamber, RS13/5 output
2045-1463		Quartz ' <b>Libra</b> ' cyclonic spray chamber for PC3 thermostatted housing
2045-1265		' <b>Electron</b> ' Mini-Cyclonic Spray Chamber, RS13/5 output
2045-1267		' <b>Positron</b> ' Baffled Mini-Cyclonic Spray Chamber, RS13/5 output
3454-1060		Quartz ' <b>Electron</b> ' Mini-Cyclonic Spray Chamber
2124-5471		<b>Coaxial</b> HF-resistant PFA spray chamber, diameter 47mm.
2029-3111	3160-3111	<b>PFA baffled</b> Cyclonic Spray Chamber, with auxiliary gas input
2029-2088		' <b>Erebus</b> ' HF-resistant Cyclonic Spray Chamber, RS13/5 output, 4mm drain
2045-0170	3203739	VG Thermo Controlled Maxim Spray Chamber
2055-0180	1040950	JC-13 Metal Spring Clamp
2045-0030	3601145	Standard Quartz Torch for X-series ICP-MS, 2 locating beads, 1.5mm Injector
2045-0040	3601147	Standard Quartz Torch for X-series ICP-MS, 2 locating beads, 2mm Injector
2045-0040	ES-1203-0010	Standard Quartz Torch for X-series ICP-MS, 2 locating beads, 1mm Injector
2045-0036	3600884	Quartz torch bonnet
2045-0250	3601149	Quartz Bonnet, X series 2
2045-0260	1200274	Quartz torch bonnet cap
2045-0038	3601149	Quartz torch bonnet
2045-0290	3203789	Torch Shield
2045-9900	3601219	PlasmaScreen Silver Screen for torch
6608-1102	EP-1202-1203	HF-res Semi-Demountable Torch with Sapphire Injector, O-ring free

EP-1203-3180		O-ring free Sapphire Injector for 6608-1102, ID = 1,8mm
EP-1202-0301		Quartz Torch for 6608-1102
EP-5509-0001		PFA spray-chamber adaptor for 6608-1102
EP-1202-0499		Complete Semi-Demountable Torch with O-ring free quartz Injector
EP-1202-0301		Precision Quartz Torch for O-ring free base
EP-1253-0150		Quartz Injector with PFA O-ring free base, ID = 1,5 mm or 1,0 mm.
2045-0055		SDT Quartz Torch for X-series ICP-MS, 2 locating beads, 1.5mm Injector
2045-0300	3600169	Quartz Elbow
2045-0320	3200727	Elbow, Ball/Socket Quartz
<b>EPOND No.</b>	<b>OEM P/No.</b>	<b>Cones – PQ, PlasmaTrace, ExCell and X-series</b>
VG1001-Ni	3004661	Nickel Sampler ( 10% discount for P/No. 4001690 - pack of 3 ** )
VG1004-Ni	3200860	Nickel Microskimmer ( 10% discount for P/No. 4200846 - pack of 3 ** )
VG1006-Pt	3200262	Platinum Sampler Cone (18mm insert)
VG1008-Pt	3201101	Platinum Micro skimmer Cone ( Enhanced Interface )
VG1012-Ni	3201585	Nickel PQE Skimmer Cone
VG5001_2	3004382	Quantity 2 of, Graphite Gasket for Sampler cones, Thermo # 1042620
<b>EPOND No.</b>	<b>OEM P/No.</b>	<b>Xi – series cones</b>
VG1021-Ni	3600812	Nickel Sampler Cone, 100% Nickel, same as VG1021-Ni-e
VG1021-Ni/Cu	3600812	Nickel Sampler Cone, Cu core / ES-3000-1221
VG1022-Ni	3600811	Nickel Skimmer Cone
VG1023-Cu	3600813	Copper Skimmer Adaptor
VG1026A-Pt/Cu	3601289	Platinum Sampler Cone with 15,2mm Pt insert
VG1026-Pt	1067501	Platinum Sampler Cone
VG1028-X-Pt	1142151	Platinum X-Skimmer Cone
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe
<b>EPOND No.</b>		<b>Accessory Description</b>
QuDIN-1011		Mermet-Todoli <b>Quasi-DIN</b> adaptor, with Peri-Pump tubes and connectors
ES-4360-1000		<b>Apex Q</b> System for XSeries2 with nebulizer, uptake capillaries, gas lines and fittings.
ES-4460-1000		<b>Apex HF</b> System for XSeries2 with nebulizer, uptake capillaries, gas lines and fittings.
SC-8106-9000		<b>SC-Micro</b> Ultra-clean autosampler with dual flowing rinse
ES-4599-3002		<b>Spiro TMD</b> Membrane Desolvation Module ( latest vsn of ACM ), with Power Supply
2019-4260		'XENA' HF-resistant micro-flow meter / aspiration monitor, range : 5 – 1999 ul/min.
2059-5260		'KEPLER' HF-resistant continuous-reading uptake meter , range : 5 – 200 ul/min.
8100-2051		'PEARL' Cleaning Device for concentric glass nebulizers
6808-2280		'ISOBAR' Internal Standard Kit, HF-resistant, with tubing and connectors
<b>EPOND N°</b>	<b>OEM P/No.</b>	<b>Description, tubes</b>
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		'CT – style' Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10
ALK-025		ArgoLock kit for 'Quick-Fit' argon connection to nebulizer
<b>ERT-040</b>		<b>EPOND 'RubeTubes', sample tube connectors, packet of 40</b>
DLK-075		DrainLink nebulizer adaptor and PTFE tube, pkt 4
NLK-030		Neb-Link sample tubes, packet of 8, L=70cm.
SAM-06A		SamTubes, torch side arm connectors, 6mm, pkt 4
		<b>Autosamplers and ICP / AAS Standard and Reference Solutions</b>
SC-8000-XXXX		<b>Complete range of ESI SC-series autosamplers and spare parts - please ask for details</b>
SC-0305-1		Autosampler 5-litre rinse bottle, for SC-DX autosamplers
LR-60-16		Rack, 60 (5x12) positions, for 16mm OD vials
MP2-4CH		Peristaltic Pump, 4-channels, autosampler peripheral
ICP-standards		
		... custom standards supplied - according to your specifications ...
		<b>www.e-pond.biz</b>