

2010		EPOND Price list
EPOND N°	OEM P/No.	Varian Vista, VARIAN 725 and Liberty ICP-OES, radial-view
2050-0710N	20-100765-00	C-Type Concentric Glass Nebulizer, 2 ml/min., with NebLink
2050-0800N	20-100695-00	K-Type Concentric Glass Nebulizer, 2 ml/min., with NebLink
2050-0710NT		EPOND 'Tyfoon' Hi-performance Concentric Glass Nebulizer, with NebLink sample tube
2050-M030-C1	TR30-C1	'Meinhard' Concentric glass nebulizer, Type C
2080-0550	99-100574-00	Varian V-groove nebulizer, 0.7 l/m Ar, 35 psi
2051-0100		Nordermeer V-groove nebulizer
2050-4250		'NB-Class' (No-Break) Boreal Concentric Polymer Nebulizer, for low uptakes
2550-2001		HF-resistant PFA Concentric Nebulizer, 80 ul/min uptake version
2550-2002		HF-resistant PFA Concentric Nebulizer, 100 ul/min uptake version
2550-2005		HF-resistant PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		PolyPro-ST HF-res concentric nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2030-2028	ES-2030-28	PFA-HS HF-res. High-Solids concentric neb, with detachable sample uptake tube
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Teflon Burgener MiraMist Parallel-Path Nebulizer
3150-0720		PEEK Burgener MiraMist Parallel-Path Nebulizer
3150-0730		Teflon T2100 'Work-horse' Parallel-Path Nebulizer
2051-0262		'Libra' Baffled Cyclonic Glass spray chamber, with polymer nebulizer adaptor
2051-0260		'Althea' Cyclonic Glass spray chamber, glass nebulizer adaptor with 2 O-rings
2051-0264		'Libra' Baffled Cyclonic Glass spray chamber, glass nebulizer adaptor with 2 O-rings
2051-0263		Jacketed 'Althea' Cyclonic Glass spray chamber, water-cooled
2080-0210E	2010081700	'Althea' Cyclonic Spray Chamber, with polymer nebulizer adaptor
2080-0212		'Electron' Mini-Cyclonic Spray Chamber 8mm Top, 4mm Drain
2080-0213		'Positron' Baffled Mini-Cyclonic Spray Chamber 8mm Top, 4mm Drain
2080-0170	9910070400	Varian Water Cooled Glass Spray Chamber
2051-026T		HF-resistant Cyclonic spray chamber, PTFE
EP-3165-3001		PFA Cyclonic Spray Chamber
EP-2128-0351	ES-2128-0351	35 mm PFA Spray Chamber and ESI 'PureCap' end cap
EP-2128-0471	ES-2128-0471	47 mm PFA Spray Chamber and ESI 'PureCap' end cap
EP-2128-1000	ES-2128-1000	Mounting Bracket for PFA Spray Chamber
EP-4157-1000	ES-4157-1000	PC3 Peltier-cooled Cyclonic Quartz Spray Chamber, with nebulizer and Power Supply
EP-9100-6600		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2080-2929		HF-resistant ' Erebus ' Cyclonic spray chamber, with drain tube
2080-0050	20-1006969-00	Varian radial-view ICPQuartz Torch, Low Flow, 1.5mm Injector
2080-005B	20-100696-00B	Varian radial-view ICP Torch with Ball Joint
2080-005W		Varian radial-view ICP Wear Metals Torch
		*** discount of 20% for orders of 3 or more torches ***
2080-0080	20-100764-00	Varian Radial Liberty Extended torch, for N2 determination

2080-0100	20-100729-00	Varian Radial Liberty, HF-res Demountable Quartz Torch Body
2080-0250	20-100707-90	Varian Radial Quartz Torch Bonnet
2080-0410	20-100730-00	Varian Radial Liberty Quartz Injector 1.8mm ID
2080-0420		Varian Radial Liberty Quartz Injector 0.8mm ID, for organics
2080-0450	24-100675-00	Varian Radial Liberty Alumina Injector 1.8mm ID
2080-0460		Varian Radial Liberty Alumina Injector 2.5mm ID
2080-0480	99-100573-00	Varian Radial Liberty Teflon Injector Holder for demountable torch
2080-0500	99-100568-00	Varian radial-view ICP-ES Demountable Torch Kit, Includes: Torch Body 2080-0100, Alumina Injector 2080-0450 & Teflon Injector Holder 2080-0480
EPOND N°		Other Accessories
APEX-Q-VAR		APEX-Q Desolvator, with PolyPro-ST nebulizer and 'Pluto' flow meter - improves sensitivity
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module (latest model of ACM), with Power Supply
2019-4260		' XENA ' micro-flow meter/monitor, range : 5 - 1999 ul/min., mains powered
8100-2051		' PEARL ' Cleaning Device for concentric glass nebulizers - suits all types concentric glass nebs
6808-2280		' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors
2045-0700		' SOLVIT ' Argon Humidifier, with tubing and connectors
EPWP-0710	WP-710	Water Purification System - produces 30 litres/hour of UP water
EPME-PG41-100		100-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm
EPME-PG41-500		500-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm
		reference solutions can be made to measure
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40
DLK-075		<i>DrainLink</i> spray chamber drain tube with 0,75mm ID, pkt 5
SAM-06A		<i>SamTubes</i> torch side arm connectors, 6mm, pkt 4.
SAM-08A		<i>SamTubes</i> torch side arm connectors, 8mm, pkt 3.
		Autosamplers and ICP / AAS Standard and Reference Solutions
SC-8000-XXXX		<i>Complete range of ESI SC-series autosamplers - please ask for details</i>
ICP-standards		Complete range of Certified Mono-Element and Multi-Element standard solutions
		*** custom standards made to your specifications ***
		www. e-pond. biz