

| 2010 | | EPOND S.A. Parts list |
|-----------------|------------------|--|
| EPOND N° | OEM P/No. | Description, parts for Varian UltraMass ICP-MS |
| 2050-0710N | 2010076500 | C-Type ' Caliber ' standard Concentric Glass Nebulizer, 2 ml/min., with 'NebLink' tube |
| 2050-0710NT | | Tyfoon High-performance concentric glass nebulizer, with NebLink tube |
| 2050-4150 | | 'NB-Class' (No-Break) Zefyr Concentric Polymer Nebulizer, low uptake version |
| 2050-4250 | | 'NB-Class' (No-Break) Boreal Concentric Polymer Nebulizer, medium uptake version |
| 2550-2000 | | HF-resistant MicroFlow PFA Concentric Nebulizer, 25 ul/min uptake version |
| 2550-2001 | | HF-resistant MicroFlow PFA Concentric Nebulizer, 80 ul/min uptake version |
| 2550-2002 | | HF-resistant MicroFlow PFA Concentric Nebulizer, 120 ul/min uptake version |
| 2550-2005 | | HF-resistant MicroFlow PFA Concentric Nebulizer, 360 ul/min uptake version |
| 2550-2028 | | PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube |
| 2550-2024 | | PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube |
| 2550-2045 | ES-2045 | Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min. |
| 2550-2043 | ES-2043 | Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min. |
| 2550-2042 | ES-2042 | Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min. |
| 2550-2041 | ES-2041 | Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min. |
| | Note : | <i>Sample capillaries above can be supplied with probes to suit all autosampler types</i> |
| 3150-0710 | | Teflon Burgener MiraMist Parallel-Path Nebulizer |
| 3150-0720 | | PEEK Burgener MiraMist Parallel-Path Nebulizer |
| | | |
| 2051-0264 | | 'Libra' Baffled Cyclonic Glass spray chamber, 8-mm vertical outlet |
| 2080-0170 | 9910070400 | Varian Water Cooled Spray Chamber |
| 2080-0210E | 2010081700 | Varian Cyclonic Spray Chamber 8mm Top, 4mm Drain |
| 2080-0212 | | Mini-Cyclonic Spray Chamber 8mm Top, 4mm Drain |
| 2080-2929 | | HF-resistant Cyclonic Spray Chamber |
| | | |
| 2080-0060 | 2010078300 | Varian Ultramass Quartz Torch 1.5mm Injector |
| 2080-0110 | 2010079300 | Varian UltraMass, Demountable Quartz Torch Body |
| 2080-0410 | 2010073000 | Varian UltraMass / Radial Liberty Quartz Injector ID=1.8mm |
| 2080-0420 | | Varian UltraMass / Radial Liberty Qtz. Injector ID=0.8mm |
| 2080-0450 | 2410067500 | Varian UltraMass / Radial Liberty Alumina Injector ID=1.8mm |
| 2080-0460 | | Varian UltraMass / Radial Liberty Alumina Injector ID=2.5mm |
| 2080-0480 | 9910057300 | Varian UltraMass / Radial Liberty Teflon Injector Holder |
| 2080-0510 | 9910071700 | Complete UltraMass Demountable Quartz Torch |
| EPOND N° | OEM P/No. | ICP-MS Detector |
| EP14566 | | UltraMass Electron Multiplier |
| | | |
| EPOND N° | | Special Accessories |
| APEX-Q | | APEX-Q Desolvator, with 2 PFA concentric nebulizers and PLUTO meter |
| 2019-4260 | | ' XENA ' micro-flow meter / aspiration monitor, range : 5 - 999 ul/min. |
| 8100-2051 | | ' PEARL ' Cleaning Device for concentric glass nebulizers |
| 6808-2280 | | ' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors |
| EP-4599-3002 | ES-4599-300A | Spiro TMD Membrane Desolvation Module (latest model of ACM), with Power Supply |
| | | |
| EPOND N° | OEM P/No. | Description, tubes |
| See web site | | All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12 |
| See web site | | Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10 |
| ERT-040 | | EPOND 'RubeTubes', sample tube connectors, packet of 40 |
| EST-030 | | PTFE Sample Tube, OD=1,3mm ID=0,30mm L=10m. |
| DLK-075 | | DrainLink nebulizer adaptor and tube with 0,75mm ID, pkt 5 |
| NLK-050 | | Neb-Link sample tubes, packet of 10, ID=0,50 mm, L=60cm. |
| NLK-030 | | Neb-Link sample tubes, packet of 10, ID=0,30 mm, L=60cm. |
| | | |
| ICP-standards | | <i>Certified ICP / AAS Standard solutions - www.e-pond.biz/standards</i> |
| | | |
| | | www. e-pond. biz |