

| 2010 | | VG PQ series, Axiom, ExCell, PlasmaTrace |
|--------------|--------------|---|
| EPOND N° | OEM P/No. | Description |
| 2050-0710N | 1040894 | EPOND ' Caliber ' standard Concentric Glass Nebulizer, with NebLink sample tube |
| 2050-0770Q | | EPOND A-type Concentric Quartz Nebulizer, with NebLink sample tube |
| 2050-M030-C1 | TR30-C1 | ' Meinhard ' Concentric glass nebulizer, Type C |
| 2050-4151 | | ' Zefyr ' 'NB-Class' (No-Break) Concentric Polymer Nebulizer, low-uptake version |
| 2050-4251 | | ' Boreal ' 'NB-Class' (No-Break) Concentric Polymer Nebulizer, medium uptake version |
| 2550-2020 | | HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 20 ul/min uptake version |
| 2550-2000 | | HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 50 ul/min uptake version |
| 2550-2002 | | HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 100 ul/min uptake version |
| 2550-2003 | | HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 200 ul/min uptake version |
| 2550-2005 | | HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 400 ul/min uptake version |
| 2550-2028 | | PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube |
| 2550-2024 | | PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube |
| 2030-2028 | ES-2030-28 | PFA-HS HF-res. High-Solids concentric neb, with detachable sample uptake tube |
| 2550-2045 | | Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min. |
| 2550-2043 | | Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min. |
| 2550-2042 | | Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min. |
| 2550-2041 | | Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min. |
| | | Capillaries with integrated probes - available for all common autosamplers |
| 2052-0300 | | PTFE Endcap for CG Nebulizer, 35/6 mm adaptor |
| | | |
| 2045-1136 | | Standard ' Electron ' 20ml. Mini-cyclonic glass spray chamber |
| 2045-1137 | | Standard ' Positron ' baffled Mini-cyclonic glass spray chamber |
| 2045-1136Q | | Quartz version of ' Electron ' Mini-cyclonic spray chamber |
| 2045-1132 | | Standard ' Althea ' 50ml. Cyclonic glass spray chamber |
| 2045-1132Q | | Quartz version of ' Althea ' 50ml. Cyclonic spray chamber |
| 2045-1142 | | Jacketed ' Althea ' 50ml. Cyclonic glass spray chamber |
| 2045-1134 | | ' Libra ' Baffled Cyclonic glass spray chamber |
| 2045-1144 | | Jacketed ' Libra ' Baffled Cyclonic glass spray chamber |
| 2045-0220 | 3204313 | VG Jacketed Spray Chamber with Impact Sphere |
| 2045-0220Q | | VG Quartz Jacketed Spray Chamber, with Impact Sphere |
| 2045-0130 | 3600170 | VG Quartz Conical Spray Chamber, for Peltier system |
| 2045-0140 | | VG Horizon Conical Spray Chamber for 6-mm OD Nebulizers |
| 2045-0160 | 3200841 | Quartz Coaxial VG Spray Chamber, Water Cooled |
| 2045-0170 | 3203739 | VG Thermo Controlled Maxim Spray Chamber |
| 2045-2929 | | HF-resistant ' Erebus ' Cyclonic Spray Chamber |
| 2045-0800 | | HF-resistant Teflon Baffled Spray Chamber |
| EP-2124-0351 | ES-2124-0351 | 35mm PFA spray chamber and PureCap End Cap for PQ2 and PQ3 |
| EP-6124-0001 | ES-6124-0001 | Spray Chamber Cooling Block, anodized Al, plus Mounting Bracket No. 6124-1002 |
| EP-6124-1002 | ES-6124-1002 | Mounting Bracket for attachment of Al cooling block to PQ2 |
| EP-6124-2002 | ES-6124-2002 | Mounting Bracket for attachment of Al cooling block to VG Axiom |
| 2055-0180 | 1040950 | JC-13 Metal Spring Clamp |
| | | |
| 2045-0050 | 3201192 | VG Standard Quartz Torch 1.5mm Injector |
| 2045-0060 | 3201732 | VG Standard Quartz Torch 2mm Injector |
| 2045-0080 | 1200080 | VG Demountable Torch Body, for 2045-0100 |
| 2045-0100 | | VG Semi-Demountable Torch, with Qtz Injector |
| 2045-0090 | | VG Quartz Injector, for 2045-0100 |
| 6608-0252 | | VG Semi-Demountable Torch with Inlet elbow and Quartz Injector |
| 6608-0244 | | Quartz Body for SDT 6608-0252 |
| 1200-0400 | | Semi-Demountable Torch with O-ring free Quartz Injector, ID = 1,5 mm |
| 1200-0600 | | HF-resistant Semi-Demountable Torch with O-ring free sapphire Injector, ID = 1,8 mm |
| 1200-0700 | | HF-resistant Semi-Demountable Torch with O-ring free platinum Injector, ID = 1,8 mm |
| 1233-0150 | | Quartz Injector with O-ring free PFA base, ID = 1,5 mm |
| 1203-0180 | | Sapphire Injector with O-ring free PFA base, ID = 1,8 mm |
| 1213-0180 | | Platinum Injector with O-ring free PFA base, ID = 1,8 mm |
| 1202-0200 | | Outer Quartz Torch for O-ring free Injectors |
| 1202-9000 | | VG Axiom Shield Torch Bonnet |
| 2045-0120 | 3203095 | VG Torch Body CR Horizon |
| 2045-0250 | | Standard VG Torch Bonnet, for PQ series, Micromass |
| 2045-0260 | 1200274 | VG Torch Bonnet Cap, Outside Diameter = 52.4mm |
| 2045-0270 | 3203025 | VG RF Insulator Shield, Outside Diameter = 140mm |
| 2045-0280 | 3204703 | VG Quartz Bonnet Axiom, O.D. = 72mm |

| | | |
|------------------|------------------|--|
| 2045-0290 | 3203789 | VG Torch Shield, for PlasmaScreen option, suits VG PQ ExCell |
| 2045-0300 | 3600169 | VG Quartz Elbow |
| 2045-0310 | | VG Elbow, Ball/Socket Pyrex |
| 2045-0320 | 3200727 | VG Elbow, Ball/Socket Quartz |
| 2045-0360 | 3007383 | VG UV Laser Sample Cell |
| 2045-0400 | 3201104 | VG Probe UV Laser Sample Cell |
| 2045-0420 | | VG Laser Sample Cell, Micro Probe |
| 2045-0430 | 3202287 | VG Quartz Sleeve |
| 2045-0440 | 3201759 | VG Quartz Tube |
| 2045-0450 | 3200176 | VG Quartz Sleeve #2 |
| 2045-0460 | 3201070 | VG Quartz Sleeve, Silica |
| 2045-0470 | 3201568 | VG Torch Tube Adapter |
| 2045-0490 | 3200004 | VG Quartz Socket Tube Adapter |
| | | |
| EPOND No. | OEM P/No. | Cones - PQ series, PlasmaTrace, ExCell |
| VG1001-Ni | 3004661 | Nickel Sampler (P/No. 4001690 for pack of 3 **) |
| VG1002-Ni | 3004866 | Nickel Mini-Skimmer Cone |
| VG1004-Ni | 3200860 | Nickel Micro-Skimmer Cone (enhanced interface) |
| | | <i>Cones above are available in aluminium, at nickel cone prices</i> |
| VG1006-Pt | 3200262 | Platinum Sampler Cone (15,2 mm Pt insert) |
| VG1006A-Pt | - | Platinum Sampler Cone (18 mm Pt insert) |
| VG1007-Pt | - | Platinum Mini-Skimmer Cone |
| VG1008-Pt | 3201101 | Platinum Micro skimmer Cone (Enhanced Interface) |
| VG1009-Pt | 3202005 | Platinum PT1 Skimmer Cone |
| VG1010-Ni | 3200456 | Nickel PT 1 Skimmer Cone |
| VG1011-Pt | 3201941 | Platinum PT2 Skimmer Cone |
| VG1012-Ni | 3201585 | Nickel PQE Skimmer Cone |
| VG1014-Ni | 3201457 | Nickel PT2 Micro Skimmer Cone |
| VG5001 | 3004382 | Graphite Gasket for Sampler cones |
| VG5001A | 3004382 | Graphite Sampler Cone Gasket - old interface with screw holes |
| EPOND No. | OEM P/No. | Xi - series cones |
| VG1021-Ni-e | - | eSPEC Nickel Sampler Cone, 100% Nickel |
| VG1021-Ni/Cu-e | 3600812 | eSPEC Nickel Sampler Cone, Cu core |
| VG1022-Ni-e | 3600811 | eSPEC Nickel Skimmer Cone |
| VG1023-Cu | 3600813 | Copper Skimmer Adaptor |
| VG1026-Pt | 1600335 | Platinum Sampler Cone |
| VG-1028-Pt | N.A. | Platinum Skimmer Cone |
| 2085-2206-05 | | Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe |
| 2085-2206-18 | | Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe |
| EPOND No. | OEM P/No. | Description of ICP-MS Detector |
| EPM14562 | | Detector for PQ1 and PQ2 (with 12-12 rods) |
| EPM14562A | | Detector for PQ1 and PQ2 (with SXP rods) |
| EPM4870V | 2027814 | Burle Channeltron detector - replacement now by Type 401A |
| EPM14562A | | Detector for PQ3, ExCell (sequential) |
| EPM14214 | AF 214 | Detector for PQ3, ExCell (simultaneous) |
| EPOND No. | | Accessory Description |
| APEX-Q-VG | | APEX-Q Desolvator, with 2 PFA concentric nebs and 'Sedna' flow monitor |
| EP-4599-3002 | ES-4599-300A | Spiro TMD Membrane Desolvation Module, with Power Supply |
| 2019-4260 | | ' XENA ' micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., AC powered |
| 8100-2051 | | ' PEARL ' Cleaning Device for concentric glass nebulizers |
| 6808-2280 | | ' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors |
| EPWP-0710 | WP-710 | Water Purification System - produces up to 30 litres/hour of UP water |
| | | reference solutions made to measure - on request |
| EPOND N° | OEM P/No. | Description, tubes |
| See web site | | All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12 |
| See web site | | Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10 |
| ERT-040 | | EPOND 'RubeTubes', sample tube connectors, packet of 40 |
| DLK-075 | | DrainLink nebulizer adaptor and tube with 0,75mm ID, pkt 5 |
| ICP-standards | | Certified ICP / AAS Standard solutions - www.e-pond.biz/standards |
| | | reference solutions made to measure - on request |
| | | www.e-pond.biz |