

# PFA Autosampler Vials and Caps

## The ultra pure solution container

Designed for the ESI SC-Series and other small volume autosamplers, these ultrapure vials are ideal for holding high purity, low volume samples, particularly where the sample matrix is an aggressive acid. The PFA cap can be press fitted onto the vials prior to analysis by the Autosampler, helping to prevent sample contamination and evaporation.



- Constructed from high purity, chemically resistant PFA-Teflon<sup>®</sup>.
- Available in range of sizes—1mL, 2mL and 4mL
- Press fit caps easily attach to the vials preventing contamination

## Ordering Information

ES-7220-0100	1mL PFA Vial (pack of 10)
ES-7220-0200	2mL PFA Vial (pack of 10)
ES-7220-0400	4mL PFA Vial (pack of 10)
ES-7220-0024	Press on cap for PFA vials (pack of 10)

Elemental Scientific Inc.  
2440 Cuming St.  
Omaha, NE 68131 USA  
[www.elementalscientific.com](http://www.elementalscientific.com)

Phone: 402.991.7800  
Fax: 402.991.7799  
Email: [esi@icpms.com](mailto:esi@icpms.com)





Elemental Scientific Inc.  
2440 Cuming St.  
Omaha, NE 68131 USA

Phone: 402.991.7800  
Fax: 402.991.7799  
Email: [esi@icpms.com](mailto:esi@icpms.com)