

PFA PureTip Pipette Tips

An essential component for ultratrace analysis.

Often overlooked as a source of sample contamination the pipette tip can come into direct contact with the sample or acids that are added to the sample.

Only PureTips can deliver

contamination free concentrated chemicals to the sample. PureTips can be reused hundreds of times and will not leach contaminants throughout their lifetime.

- Moulded from high purity PFA-Teflon[®], minimizing contamination
- Chemically resistant to virtually all sample types, ideal for strong acids, alkalis, organics and biologicals.
- Long life, re-usable.
- Use with standard Eppendorf type 1000 μ L & 100 μ L pipettors.

Ordering Information

ES-7000-0101 100 μ L pre-cleaned tips (pack of 10)

ES-7000-1001 1000 μ L pre-cleaned tips (pack of 10)



1000 μ L & 100 μ L moulded PFA PureTip pipette tips.

Elemental Scientific Inc.
2440 Cuming St.
Omaha, NE 68131 USA
www.elementalscientific.com

Phone: 402.991.7800
Fax: 402.991.7799
Email: esi@icpms.com

