

New PFA Vials for small volume samples

Designed for the ESI SC-Series autosamplers and other small volume autosamplers, these ultrapure vials are ideal for holding high purity, low volume samples.

- Constructed from high purity, chemically resistant PFA fluoropolymer
- Available in range of sizes
 - 1.5mL
 - 2.5mL
 - 4 mL
- Graduated volume markings
- Fin design at base allows vials to be removed easily from acid cleaning baths using forceps or tweezers
- Translucent at bottom for easy observation of sample remaining
- Conical bottom for complete sampling of micro samples
- One rack for all vial sizes
- Optional press fit caps for preventing evaporation or contamination



Specifications and Ordering Information:

Part Number	Volume	Height
V-14-0311-CG	1.5 mL	24mm
V-14-0312-CG	2.5 mL	36mm
V-14-0313-CG	4 mL	52mm

All these vials have 14mm outer diameter and ESI Micro racks for SC-autosamplers:



MR-21-14 - 21 position 3x7 rack
MR-40-14 - 4x10 rack (not shown)