

# SC-2 CLEAN ENCLOSURE

The Clean Enclosure has been developed for the SC-2 ultraclean autosampler. This enclosure helps to maintain sample integrity while samples await analysis.

- ✓ Protects samples from airborne contamination.
- ✓ Clear acrylic construction, provides sturdy protection to the SC-2.
- ✓ Covers only the sample platform
- ✓ Large door opening for easy access to the sample racks.
- ✓ Clean Enclosure door easily opens and locks away closed for rapid access to samples.
- ✓ ULPA filtration available

P/N ES-1207 Sample Enclosure  
for SC-2

Elemental Scientific Inc.  
2440 Cuming St.  
Omaha, NE 68131 USA  
[www.elementalscientific.com](http://www.elementalscientific.com)

Phone: 402.991.7800  
Fax: 402.991.7799  
Email: [esi@icpms.com](mailto:esi@icpms.com)



SC-2 autosampler with clean enclosure attached



Enclosure doors easily open for access to sample racks

