

Options

- **Clean Enclosure**

- Protects samples, equipment and personnel. Available with or without ULPA filtration.

- **FAST- Sample injection system**

- Increase sample throughput up to 5 times


- **Rinse pump**

- Teflon Diaphragm pump
- Pressurized/Gravity with polypropylene bottles or all PFA rinse path for ultra clean applications (includes PFA bottles)
- MP² – Integrated precision micro peristaltic pump



SC-E2 DXS Enclosure on Thermo ELEMENT 2

ULPA filter side plate available upon request for either the right or left side of the SC-E2 DXS enclosure



Part Number	With ULPA Filter	Description
SC-1307	SC-1307-1000	Enclosure for SC-E2 DXS Autosampler

SC-E2 DXS Autosampler

- Includes:** 2 Sample Probes
 Standards Rack (10 positions)
 50 mL Standards Vials (10 ea)
 4 Large Racks (21, 40, 60, 90 positions)
 15 mL Sample Vials (60 ea)

Part Number	Description
SC-8103-1000	SC-E2 DXS with dual flowing rinse


SC-E2 DXS FAST Autosampler

- Includes:** FAST Valve
 PFA-ST MicroFlow Nebulizer
 Standards Rack (10 positions)
 4 Large Racks (21, 40, 60, 90 positions)
 2 Sample Probes
 50 mL Standards Vials (10 ea)
 15 mL Sample Vials (60 ea)

PC³ Peltier cooled cyclonic chamber, demountable quartz torch and sapphire injector (*Complete System Only*)

Part Number	Description
E2F-52C	SC-E2 DXS FAST PC ³ Complete for ELEMENT 2
E2F-A	SC-E2 DXS FAST Basic system for ELEMENT 2 (without PC ³) Uses ELEMENT 2 spray chamber and injector

Available Sample Racks



Large Racks



Super Racks



Micro Racks



Microtiter Plates

Dimensions of the SC-E2 DXS Autosampler on the Element 2 (Optional ULPA Filter)

