

Apex Shelves for PE Elan


Greater work area for sample introduction systems

Increase the work area and make the use of sample introduction systems more convenient with these add-on shelves for PerkinElmer Elan ICPMS. Ideal for use with the Apex inlet system.

- ▲ Positions the Apex system to align with ICP torch
- ▲ Attaches to the Elan using manufacturer's contact points
- ▲ Quick and easy to install
- ▲ Provides sturdy work area capable of supporting heavy loads
- ▲ Constructed from chemically resistant polypropylene



ES-2999-5421
Shelf for PE Elan with cassette torch



ES-2999-5431
Shelf for PE Elan with bayonet torch